# Duplex sequencing for mutagenesis testing in wild-type rodents and a common human cell line Devon M. Fitzgerald<sup>1</sup>, Thao Huynh<sup>1</sup>, Raul Burciaga<sup>1</sup>, Camila Zanette<sup>1</sup>, Jake Higgins<sup>1</sup>, and Jesse Salk<sup>1</sup>

<sup>1</sup>TwinStrand Biosciences, Inc., Seattle, WA





# Introduction

Toxic compounds can alter the genetic makeup of organisms, generating mutations that cause cancer and other poor health outcomes. Standard mutagenesis assays yield a limited amount of indirect information about a compound's mutagenic potential and/or are complex and cumbersome.

**Duplex sequencing (DS)** is an error-corrected sequencing method that can directly detect, quantify, and characterize mutagenesis.

Here, we use TwinStrand DuplexSeq<sup>™</sup> Mutagenesis Assays to characterize the mutagenicity of benzo[a]pyrene (B[a]P) and N-ethyl-N-nitrosourea (ENU) across multiple tissues of wild-type rodents and in a common cell line. Also, we generated abundant, characterized source material (genomic DNA) that can be shared with the community for inter-lab reproducibility studies.

#### Mutagen exposure protocols



Male **Sprague Dawley rats** (n=3 per group) and C57BL/6Ntac mice (n=5 per group) were administered **B[a]P** or **ENU** by oral gavage, according to the schedule illustrated to the left. Untreated rodents were maintained as controls. Rodents were observed daily and kidney, liver, lung, and spleen tissue were harvested at necropsy on day 31.

#### Human TK6 lymphoblastoid cells (n= 5 dishes per treatment) were exposed to **B[a]P + S9** or **ENU** as

# **Tissue-specific mutation frequencies (rat)**



This study paves the way for streamlined genotoxicology testing, reducing animal usage, cost, and time to insight.



shown (left), or exposed to vehicle only. Doses were selected based on a range-finding study which targeted 30-60% cytotoxicity (CellTiter-Glo viability assay) and allowed for adequate growth during the post-exposure period.

**TwinStrand DuplexSeq<sup>™</sup> Mutagenesis Assays** 



**Mutagenesis Panels:** • 20 x 2.4 kb regions throughout genome • genic and intergenic regions included • regions not predicted to be involved in positive or negative selection • available for **rat** (shown left), **mouse**, and **human** 

control

B[a]P ENU

Duplex Sequencing (DS) libraries were prepared using 500 ng genomic DNA extracted from rodent tissues or cultured cells. The Enzymatic Fragmentation Module (currently in early release) was used to construct all libraries, unless otherwise noted. Hybrid capture was performed with species-matched Mutagenesis Panels.

# **Duplex sequencing performance metrics**



Both B[a]P and ENU treatments increase MF significantly across all rat tissues. The magnitude of mutagenic effects varies between tissues, with B[a]P causing the most mutagenesis in the liver and ENU having the most dramatic effect in spleen. \*\*\*  $p < 2 \times 10^{-12}$ , logistic regression; error bars = 95% CI.

# Single base-pair substitution spectra (rat)



Grouping samples by treatment and tissue reveals differences in single base-pair substitution spectra (left). treatment increases the B[a]P of C>A mutations. proportion **ENU**-treated samples are dominated Subtype by T>A and T>C mutations. Spectra C>A are consistent with previous reports.

C>G C>T Unsupervised clustering perfectly groups samples by treatment group in T>A T>C rat kidney (below left) and liver (not T>G shown) and clearly separates ENU from other treatments in lung and spleen (below right and not shown).

> Mouse and cell line datasets show similar single base-pair substitution spectra, with perfect clustering by treament (data not shown).



#### **Sequencing Errors Obscure Truth**





**Next-Generation** Sequencing (NGS)



Duplex Sequencing

ΝΊΝΊΝΙΝΝΝΊΝΙ

NNNNNNN

# **TwinStrand Duplex Sequencing<sup>™</sup> Technology** DuplexSeq<sup>™</sup> Tag

A DuplexSeq<sup>™</sup> Adapter has:

- Identical (or relatable) degenerate tags in each strand.
- An asymmetry allowing independent strand identification.



# Measuring genotoxicity



Tumor formation is triggered by exposure to





## **Enzymatic fragmentation "brightens" signal**



DS libraries were prepared using genomic DNA fragmented mechanically (Covaris) or with the new Enzymatic Fragmentation Module (early release). Enzymatic fragmentation reduces MFs of control samples, improving separation between control and treated groups and "brightening" the signal of mutagen exposure. Each data point represents a single tissue from an individual animal.

#### **Trinucleotide spectra (rat)**



ENU DS yields sufficient data to generate trinucleotide T>C T>G T>A spectra. These unique "fingerprints" of mutagenesis consider the nucleotides flanking a mutated base and can provide clues to the underlying mechanisms of mutagenesis. These "fingerprints" can be compared to databases of trinucleotide signatures identified in tumors, such as COSMIC. Trinucelotide spectra shown here represent data combined across individuals and tissues. **Trinucleotide Context** 

**Abbreviations** 

**DS** = Duplex sequencing **MF** = mutation frequency

**B[a]P** = benzo[a]pyrene **ENU** = N-ethyl-N-nitrosourea

For Research Use Only. Not for use in diagnostic procedures. ©2022 TwinStrand Biosciences, Inc. All rights reserved. All trademarks are the property of TwinStrand Biosciences, Inc. or their respective owners.

#### Mutation frequencies by species and treatment



DS reveals clear differences in MF between treatment groups in all sample types. Absolute MFs are similar between rat and mouse. TK6 cells show a 10-fold higher basal MF than rodent tissues, but this does not obscure the mutagenic effects of B[a]P or ENU. Each data point represents a single tissue from an individual animal and values are plotted on a log scale. Error bars = 95% confidence intervals (CI).

### Conclusions

 DuplexSeq<sup>™</sup> Mutagenesis Assays yield robust mutagenicity data from wild-type rodents and a common cell line after short-term exposure to two mutagens.

baseline mutation frequencies < 1-in-10-million among rodent controls</li>

 can detect subtle changes in mutation frequency (~2-fold change) and spectrum • single-base and trinucleotide spectra provide mechanistic insight into mutational processes

 demonstrates potential to use same cohort of wild-type rodents for mutagenicity testing and other routine safety testing, reducing overall animal usage

• New Enzymatic Fragmentation Module boosts assay performance by "brightening" the signal of subtle mutagen exposure.

• Source material (genomic DNA) from this study is abundant and can be shared with the community for use in inter-lab validation studies.

Contact: https://twinstrandbio.com/contact/ dfitzgerald@twinstrandbio.com